

## Product datasheet

# Human ADAM8 peptide ab22912

### 製品の概要

---

製品名 Human ADAM8 peptide

### 製品の詳細

---

由来 Synthetic

#### アミノ酸配列

生物種 Human

配列 C-QRKQGAGAPTAP

領域 813 to 824

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab22912** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-ADAM8 antibody ([ab2475](#))

製品の状態 Liquid

### 前処理および保存

---

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### 関連情報

---

機能 Possible involvement in extravasation of leukocytes.

組織特異性 Expressed on neutrophils and monocytes.

配列類似性 Contains 1 disintegrin domain.  
Contains 1 EGF-like domain.  
Contains 1 peptidase M12B domain.

細胞内局在 Membrane.

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors