

## Product datasheet

# Human DDAH1 peptide ab99047

### 1 References

#### 製品の概要

**製品名** Human DDAH1 peptide

#### 製品の詳細

**由来** Synthetic

#### アミノ酸配列

**生物種** Human

**配列** TCCSVLINKKVDS

**領域** 272 to 284

#### 特性

Our [Abpromise guarantee](#) covers the use of **ab99047** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション** Blocking - Blocking peptide for Anti-DDAH1 antibody ([ab2231](#))

**製品の状態** Lyophilised

#### 前処理および保存

**保存方法および安定性** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

**再構成** Reconstitute with water

#### 関連情報

**機能** Hydrolyzes N(G),N(G)-dimethyl-L-arginine (ADMA) and N(G)-monomethyl-L-arginine (MMA) which act as inhibitors of NOS. Has therefore a role in the regulation of nitric oxide generation.

**組織特異性** Detected in brain, liver, kidney and pancreas, and at low levels in skeletal muscle.

**配列類似性** Belongs to the DDAH family.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors