

Product datasheet

Human PACSIN1 peptide ab45457

製品の概要

製品名 Human PACSIN1 peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 C-SSSYDEASLAPEET

領域 2 to 15

特性

Our [Abpromise guarantee](#) covers the use of **ab45457** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-PACSIN1 antibody ([ab2532](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

関連情報

機能 May play a role in vesicle formation and transport.

組織特異性 Highly expressed in brain and in a lesser extent in heart, pancreas and liver.

配列類似性 Belongs to the PACSIN family.
Contains 1 FCH domain.
Contains 1 SH3 domain.

翻訳後修飾 Phosphorylated by casein kinase 2 (CK2) and protein kinase C (PKC).

細胞内局在 Cytoplasm. Colocalizes with dynamin I at vesicular structures in the cell body and neurites.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors