

## Product datasheet

# Human CASP peptide ab22893

### 製品の概要

**製品名** Human CASP peptide

### 製品の詳細

**由来** Synthetic

**アミノ酸配列**

**生物種** Human

**配列** SLQRLQHSNGN-C

**領域** 2 to 14

### 特性

Our [Abpromise guarantee](#) covers the use of **ab22893** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション** Blocking - Blocking peptide for Anti-CASP antibody ([ab2247](#))

**製品の状態** Liquid

**備考** Protein previously labeled as PSCDBP.

### 前処理および保存

**保存方法および安定性** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### 関連情報

**機能** By its binding to cytohesin-1 (CYTH1), it modifies activation of ARFs by CYTH1 and its precise function may be to sequester CYTH1 in the cytoplasm.

**組織特異性** Expressed in lymph nodes, thymus, spleen, lung, peripheral blood leukocytes and bone marrow.

**配列類似性** Contains 1 PDZ (DHR) domain.

**細胞内局在** Cytoplasm. Early endosome. Recruited from the cytosol to endosomes by SNX27.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors