

## Product datasheet

# Human Syntrophin gamma 2 peptide ab23267

### 製品の概要

---

**製品名** Human Syntrophin gamma 2 peptide

### 製品の詳細

---

**由来** Synthetic

#### アミノ酸配列

**生物種** Human

**配列** C-DSQSLARKYMYSS

**領域** 527 to 539

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab23267** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション** Blocking - Blocking peptide for Anti-Syntrophin gamma 2 antibody ([ab15713](#))

**製品の状態** Liquid

### 前処理および保存

---

**保存方法および安定性** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### 関連情報

---

**機能** Adapter protein that binds to and probably organizes the subcellular localization of a variety of proteins. May link various receptors to the actin cytoskeleton and the dystrophin glycoprotein complex.

**組織特異性** Widely expressed. Strong expression in brain and testis. In CNS, it is expressed in the perikaryon and proximal portion of the neuronal processes. Strong expression in the hippocampus, neuron-rich dentate granule cells, and pyramidal cell layers. Highly expressed in neurons of the cerebral cortex. Also expressed in the cerebellar cortex, deep cerebellar nuclei, thalamus, and basal ganglia.

配列類似性	Belongs to the syntrophin family. Contains 1 PDZ (DHR) domain. Contains 1 PH domain.
ドメイン	The association with dystrophin or related proteins probably leaves the PDZ domain available to recruit proteins to the membrane.
細胞内局在	Cell membrane > sarcolemma. Cytoplasm > cytoskeleton. In skeletal muscle, it localizes at the cytoplasmic side of the sarcolemmal membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors