

Product datasheet

Human TIGD1 peptide ab23179

製品の概要

製品名 Human TIGD1 peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 C-PAKRVRLTEGSD

領域 580 to 591

特性

Our [Abpromise guarantee](#) covers the use of **ab23179** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-TIGD1 antibody ([ab10092](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

関連情報

配列類似性 Belongs to the tigger transposable element derived protein family.
Contains 1 DDE domain.
Contains 1 HTH CENPB-type DNA-binding domain.
Contains 1 HTH psq-type DNA-binding domain.

細胞内局在 Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors