

## Product datasheet

# Human gamma H2A.X (phospho S139) peptide ab15645

1 References [画像数 1](#)

### 製品の概要

**製品名** Human gamma H2A.X (phospho S139) peptide

### 製品の詳細

**由来** Synthetic

#### アミノ酸配列

**アクセッション番号** [P16104](#)

**生物種** Human

**修飾** phospho S139

### 特性

Our [Abpromise guarantee](#) covers the use of **ab15645** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション** Blocking - Blocking peptide for Anti-gamma H2A.X (phospho S139) antibody - ChIP Grade ([ab2893](#))

**製品の状態** Liquid

**備考**

- First try to dissolve a small amount of peptide in either water or buffer. The more charged residues on a peptide, the more soluble it is in aqueous solutions.
- If the peptide doesn't dissolve try an organic solvent e.g. DMSO, then dilute using water or buffer.
- Consider that any solvent used must be compatible with your assay. If a peptide does not dissolve and you need to recover it, lyophilise to remove the solvent.
- Gentle warming and sonication can effectively aid peptide solubilisation. If the solution is cloudy or has gelled the peptide may be in suspension rather than solubilised.
- Peptides containing cysteine are easily oxidised, so should be prepared in solution just prior to use.

### 前処理および保存

**保存方法および安定性** Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

## 関連情報

---

機能	<p>Variant histone H2A which replaces conventional H2A in a subset of nucleosomes. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling. Required for checkpoint-mediated arrest of cell cycle progression in response to low doses of ionizing radiation and for efficient repair of DNA double strand breaks (DSBs) specifically when modified by C-terminal phosphorylation.</p>
配列類似性	<p>Belongs to the histone H2A family.</p>
発生段階	<p>Synthesized in G1 as well as in S-phase.</p>
ドメイン	<p>The [ST]-Q motif constitutes a recognition sequence for kinases from the PI3/P4-kinase family.</p>
翻訳後修飾	<p>Phosphorylated on Ser-140 (to form gamma-H2AFX or H2AX139ph) in response to DNA double strand breaks (DSBs) generated by exogenous genotoxic agents and by stalled replication forks, and may also occur during meiotic recombination events and immunoglobulin class switching in lymphocytes. Phosphorylation can extend up to several thousand nucleosomes from the actual site of the DSB and may mark the surrounding chromatin for recruitment of proteins required for DNA damage signaling and repair. Widespread phosphorylation may also serve to amplify the damage signal or aid repair of persistent lesions. Phosphorylation of Ser-140 (H2AX139ph) in response to ionizing radiation is mediated by both ATM and PRKDC while defects in DNA replication induce Ser-140 phosphorylation (H2AX139ph) subsequent to activation of ATR and PRKDC. Dephosphorylation of Ser-140 by PP2A is required for DNA DSB repair. In meiosis, Ser-140 phosphorylation (H2AX139ph) may occur at synaptonemal complexes during leptotene as an ATM-dependent response to the formation of programmed DSBs by SPO11. Ser-140 phosphorylation (H2AX139ph) may subsequently occurs at unsynapsed regions of both autosomes and the XY bivalent during zygotene, downstream of ATR and BRCA1 activation. Ser-140 phosphorylation (H2AX139ph) may also be required for transcriptional repression of unsynapsed chromatin and meiotic sex chromosome inactivation (MSCI), whereby the X and Y chromosomes condense in pachytene to form the heterochromatic XY-body. During immunoglobulin class switch recombination in lymphocytes, Ser-140 phosphorylation (H2AX139ph) may occur at sites of DNA-recombination subsequent to activation of the activation-induced cytidine deaminase AICDA. Phosphorylation at Tyr-143 (H2AXY142ph) by BAZ1B/WSTF determines the relative recruitment of either DNA repair or pro-apoptotic factors. Phosphorylation at Tyr-143 (H2AXY142ph) favors the recruitment of APBB1/FE65 and pro-apoptosis factors such as MAPK8/JNK1, triggering apoptosis. In contrast, dephosphorylation of Tyr-143 by EYA proteins (EYA1, EYA2, EYA3 or EYA4) favors the recruitment of MDC1-containing DNA repair complexes to the tail of phosphorylated Ser-140 (H2AX139ph).</p> <p>Monoubiquitination of Lys-120 (H2AXK119ub) by RING1 and RNF2/RING2 complex gives a specific tag for epigenetic transcriptional repression. Following DNA double-strand breaks (DSBs), it is ubiquitinated through 'Lys-63' linkage of ubiquitin moieties by the E2 ligase UBE2N and the E3 ligases RNF8 and RNF168, leading to the recruitment of repair proteins to sites of DNA damage. Monoubiquitination and ionizing radiation-induced 'Lys-63'-linked ubiquitination are distinct events.</p>
細胞内局在	<p>Nucleus. Chromosome.</p>



Rabbit polyclonal to H2A.X (phospho S139) - [ab2893](#) used at 1/500.

Blocking peptides at 1ug/lane.

Lanes 1, 3 and 5: Control HeLa cell histone preparation

Lanes 2, 4 and 6: Colcemid treated HeLa histone preparation

Lanes 1 and 2: [ab2893](#)

Lanes 3 and 4: [ab2893](#) + Histone H2A.X (phospho Ser 139) peptide - [ab15645](#)

Lanes 5 and 6: [ab2893](#) + Histone H2A.X non-phospho peptide - [ab15646](#)

Secondary antibody - [ab6721](#) - Goat polyclonal to Rabbit IgG (HRP) at 1/5000.

Rabbit polyclonal to H2A.X (phospho S139) - [ab2893](#) used at 1/500.

Blocking peptides at 1ug/lane.

Lanes 1, 3 and 5: Control HeLa cell histone preparation

Lanes 2, 4 and 6: Colcemid treated HeLa histone preparation

Lanes 1 and 2: [ab2893](#)

Lanes 3 and 4: [ab2893](#) + Histone H2A.X (phospho Ser 139)peptide - [ab15645](#)

Lanes 5 and 6: [ab2893](#) + Histone H2A.X non-phospho peptide - [ab15646](#)

Secondary antibody - [ab6721](#) - Goa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors