


Product datasheet

Anti-ICAM4 antibody ab112553

画像数 1

製品の概要

製品名	Anti-ICAM4 antibody
製品の詳細	Rabbit polyclonal to ICAM4
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to amino acids 43-210 of Human ICAM4 (BC000046).
ポジティブ・コントロール	Human fetal colon and fetal brain lysates.

法規制情報

医薬用外毒物

製品の特性

製品の状態	Lyophilised: Reconstitute with 200ul distilled sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituent: 99% PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112553** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 29 kDa (predicted molecular weight: 29 kDa).

ターゲット情報

機能	ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). ICAM4 is also a ligand for alpha-4/beta-1 and alpha-V integrins.
組織特異性	Erythrocytes.
配列類似性	Belongs to the immunoglobulin superfamily. ICAM family. Contains 2 Ig-like C2-type (immunoglobulin-like) domains.
翻訳後修飾	N- and O-glycosylated.
細胞内局在	Secreted and Cell membrane.

画像



All lanes : Anti-ICAM4 antibody (ab112553)
at 1/500 dilution

Lane 1 : Human fetal colon lysate

Lane 2 : Human fetal brain lysate

Predicted band size: 29 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors