

Product datasheet

Anti-IL-17D antibody ab22063

画像数 1

製品の概要

製品名	Anti-IL-17D antibody
製品の詳細	Rabbit polyclonal to IL-17D
由来種	Rabbit
特異性	ab22063 reacts with IL-27/ IL-30.
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Human
免疫原	Synthetic peptide: L LVQAGVWGFP RPPGRPQLSL QELRREFTVS LHLARKLLSE VRGQHRFAE , corresponding to N terminal amino acids 20-70 of Human IL27. Run BLAST with Run BLAST with
ポジティブ・コントロール	Jurkat, NIH 3T3

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium azide Constituents: PBS, 0.05% BSA
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab22063** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.001 - 0.002 mg/ml. Detects a band of approximately 27 kDa.

ターゲット情報

機能	Induces expression of IL-6, IL-8, and GM-CSF from endothelial cells.
組織特異性	Expressed preferentially in adipose, skeletal muscle and CNS.
配列類似性	Belongs to the IL-17 family.
細胞内局在	Secreted.

画像



All lanes : Anti-IL-17D antibody (ab22063)

Lane 1 : Identification of IL-27 in Jurkat cells

Lane 2 : Identification of IL-27 in NIH-3T3 cells

Lysates/proteins at 1 µg/ml per lane.

Observed band size: 27 kDa

[why is the actual band size different from the predicted?](#)

Additional bands at: 45 kDa, 50 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors