

Product datasheet

Anti-IL1RA antibody ab60978

画像数 1

製品の概要

製品名	Anti-IL1RA antibody
製品の詳細	Sheep polyclonal to IL1RA
由来種	Sheep
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant full length protein (Human)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.0536% PBS, 0.87% Sodium chloride
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab60978** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

追加情報

WB: Use at an assay dependent dilution. Detects bands of approximately 20 and 40 kDa (predicted molecular weight: 20 kDa). We have no information regarding the 40 kDa band.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	Inhibits the activity of interleukin-1 by binding to receptor IL1R1 and preventing its association with the coreceptor IL1RAP for signaling. Has no interleukin-1 like activity. Binds functional interleukin-1 receptor IL1R1 with greater affinity than decoy receptor IL1R2; however, the physiological relevance of the latter association is unsure.
組織特異性	The intracellular form of IL1RN is predominantly expressed in epithelial cells.
関連疾患	Microvascular complications of diabetes 4 Interleukin 1 receptor antagonist deficiency
細胞内局在	Cytoplasm and Secreted.

画像



Anti-IL1RA antibody (ab60978) at 5 µg/ml +
Recombinant human IL-1Ra. at 0.1 µg

Secondary

anti-sheep IgG - HRP at 1/10000 dilution

Predicted band size: 20 kDa

Observed band size: 20 kDa

Additional bands at: ~40 kDa (possible dimer)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors