

Product datasheet

Anti-IL-3RA / CD123 antibody [6H6] (Allophycocyanin) ab95781

画像数 1

製品の概要

製品名	Anti-IL-3RA / CD123 antibody [6H6] (Allophycocyanin)
製品の詳細	Mouse monoclonal [6H6] to IL-3RA / CD123 (Allophycocyanin)
由来種	Mouse
標識	Allophycocyanin. Ex: 645nm, Em: 660nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to Human IL-3RA/ CD123.
ポジティブ・コントロール	Normal Human peripheral blood cells
特記事項	Previously labelled as IL3RA.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS, BSA, Gelatin, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
クローン名	6H6
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab95781** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 5µl for 10 ⁶ cells.

ターゲット情報

機能	This is a receptor for interleukin-3.
配列類似性	Belongs to the type I cytokine receptor family. Type 5 subfamily.
ドメイン	The WSXWS motif appears to be necessary for proper protein folding and thereby efficient intracellular transport and cell-surface receptor binding. The box 1 motif is required for JAK interaction and/or activation.
細胞内局在	Membrane.

画像



Staining of normal human peripheral blood cells with mouse IgG1 k Isotype Control APC (open histogram) or ab95781 (filled histogram). Cells in the lymphocyte gate were used for analysis.

Flow Cytometry - Anti-IL-3RA / CD123 antibody
[6H6] (Allophycocyanin) (ab95781)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors