

Product datasheet

Anti-Interferon alpha 2b antibody [24B10] ab87900

画像数 1

製品の概要

製品名	Anti-Interferon alpha 2b antibody [24B10]
製品の詳細	Mouse monoclonal [24B10] to Interferon alpha 2b
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide derived Human Interferon alpha 2b protein. Immunogen の所有権に関して
ポジティブ・コントロール	This antibody gave a positive signal against a Interferon alpha 2b recombinant protein.
特記事項	<p>This antibody clone is manufactured by Abcam.</p> <p>If you require this antibody in a particular buffer formulation or a particular conjugate for your experiments, please contact orders@abcam.com or you can find further information here.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 1% BSA, PBS. pH 7.4
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	24B10
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab87900** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 5 µg/ml. Detects a band of approximately 21 kDa (predicted molecular weight: 21 kDa).

ターゲット情報

機能	Produced by macrophages, IFN-alpha have antiviral activities.
配列類似性	Belongs to the alpha/beta interferon family.
細胞内局在	Secreted.

画像



Anti-Interferon alpha 2b antibody [24B10] (ab87900) at 5 µg/ml + HeLa expressing Interferon alpha 2b recombinant protein. at 10 µg

Secondary

Goat Anti-Mouse IgG H&L (HRP) preadsorbed (ab97040) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 21 kDa

Observed band size: 21 kDa

Exposure time: 3 minutes

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors