

## Product datasheet

# Anti-L Ficolin antibody [16] ab112558

### 製品の概要

製品名	Anti-L Ficolin antibody [16]
製品の詳細	Mouse monoclonal [16] to L Ficolin
由来種	Mouse
特異性	No cross reactivity is seen with M Ficolin (Ficolin 1) and H Ficolin (Ficolin 3).
アプリケーション	<b>適用あり:</b> Sandwich ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Purified procaryote expressed recombinant Human L Ficolin and full-length L Ficolin (Ficolin 2) protein adsorbed onto aluminum hydroxide gel.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 2.9% Sodium chloride, 97% Phosphate Buffer
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	16
ミエローマ	Sp2
アイソタイプ	IgG1
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112558** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Sandwich ELISA		1/2000. as capture antibody.
WB		Use at an assay dependent concentration. Predicted molecular weight: 34 kDa.

## ターゲット情報

**関連性** May function in innate immunity through activation of the lectin complement pathway. Calcium-dependent and GlcNAc-binding lectin. Enhances phagocytosis of *S.typhimurium* by neutrophils, suggesting an opsonic effect via the collagen region.

**細胞内局在** Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors