

## Product datasheet

# Anti-LRRC56 antibody ab122506

### 画像数 2

#### 製品の概要

製品名	Anti-LRRC56 antibody
製品の詳細	Rabbit polyclonal to LRRC56
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, ICC/IF
種交差性	<b>交差種:</b> Human
免疫原	antigen, corresponding to amino acids 29-133 of Human LRRC56 (Q8Y6G6).
ポジティブ・コントロール	Human testis tissue.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122506** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.


アプリケーション	Abreviews	特記事項
----------	-----------	------

ICC/IF                      Use a concentration of 1 - 4 µg/ml.  
Recommend PFA Fixation and Triton X-100 treatment

### ターゲット情報


**配列類似性**                      Contains 5 LRR (leucine-rich) repeats.  
Contains 1 LRRCT domain.

### 画像



Immunofluorescent staining of Human cell line U-251 MG shows positivity in nucleus but not nucleoli. Recommended concentration of ab122506 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-LRRC56 antibody (ab122506)



ab122506, at 1/200 dilution, staining LRRC56 in formalin-fixed, paraffin-embedded Human testis tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-LRRC56 antibody (ab122506)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors