

Product datasheet

Anti-LRRC67 antibody ab122006

画像数 3

製品の概要

製品名	Anti-LRRC67 antibody
製品の詳細	Rabbit polyclonal to LRRC67
由来種	Rabbit
アプリケーション	適用あり: WB, ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	antigen sequence corresponding to amino acids 1-73 (MVRLTLDLIA RNSNLKPRKE ETISQCLKKI THINFSDKNI DAIEDLSLCK NLSVLYLYDN CISQITNLNY ATN) of Human LRRC67. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human stomach tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122006** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/250. Predicted molecular weight: 41 kDa.


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment
IHC-P		1/800. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

配列類似性


Contains 7 LRR (leucine-rich) repeats.
Contains 1 LRRCT domain.

画像



Immunofluorescent staining of Human cell line U-2 OS shows positivity in nucleus but not nucleoli and cytoplasm. Recommended concentration of ab122006 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-LRRC67 antibody (ab122006)



Western blot - Anti-LRRC67 antibody (ab122006)

All lanes : Anti-LRRC67 antibody (ab122006)
at 1/100 dilution

Lane 1 : RT-4 cell lysate


Lane 2 : U-251 MG cell lysate

Lane 3 : Human plasma

Lane 4 : Human liver

Lane 5 : Human tonsil

Predicted band size: 41 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-LRRC67 antibody (ab122006)

ab122006, at a 1/800 dilution, staining LRRC67 in paraffin-embedded Human stomach tissue by Immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors