

## Product datasheet

# Mouse monoclonal [SB62a] Anti-Rabbit IgG light chain ab99695

### 製品の概要

製品名	Mouse monoclonal [SB62a] Anti-Rabbit IgG light chain
由来種	Mouse
ターゲット生物種	Rabbit
アプリケーション	<b>適用あり:</b> WB, ELISA
免疫原	Rabbit Light Chain

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 8.20 Constituent: 100% Borate buffered saline
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	SB62a
アイソタイプ	IgG2b
軽鎖の種類	lambda

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99695** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.
ELISA		Use a concentration of 1 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors