

Product datasheet

Anti-NYS48 antibody ab95970

画像数 4

製品の概要

製品名	Anti-NYS48 antibody
製品の詳細	Rabbit polyclonal to NYS48
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P, ICC/IF
種交差性	交差種: Human
免疫原	The immunogen is a recombinant protein corresponding to a region within amino acids 216-410 of human NYS48
ポジティブ・コントロール	K562 and HeLa lysates

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 1.21% Tris, 0.75% Glycine
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab95970** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/2000. Predicted molecular weight: 45 kDa.
IHC-P		1/50 - 1/500.
ICC/IF		1/100 - 1/200.

ターゲット情報

画像




All lanes : Anti-NYS48 antibody (ab95970) at 1/500 dilution


Lane 1 : His-Hice 1 (2x)

Lane 2 : His-Hice 1 (1x)

Predicted band size: 45 kDa




Immunohistochemical analysis of paraffin-embedded human lymph tissue using ab95970 at 1/50 dilution. (Image 20x)



Immunohistochemical analysis of paraffin-embedded human lymph tissue using ab95970 at 1/50 dilution. (Image 40x)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NYS48 antibody (ab95970)



Immunofluorescence analysis of human osteosarcoma cell line U2OS using ab95970 against NYS48 (1), alpha tubulin (2) and merged with DAPI (3) at 1/50-1/200 dilution

Immunocytochemistry/ Immunofluorescence - Anti-NYS48 antibody (ab95970)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors