

## Product datasheet

# Anti-PCNP antibody [2707C2a] ab70733

### 画像数 1

#### 製品の概要

製品名	Anti-PCNP antibody [2707C2a]
製品の詳細	Mouse monoclonal [2707C2a] to PCNP
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, Dot blot
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment derived from internal sequence of human PCNP
ポジティブ・コントロール	Recombinant human PCNP.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
精製度	Protein G purified
特記事項 (精製)	ab70733 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum and filtered through a 0.22µm membrane.
ポリ/モノ	モノクローナル
クローン名	2707C2a
アイソタイプ	IgG1

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab70733** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 19 kDa.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.

## ターゲット情報

機能	May be involved in cell cycle regulation.
翻訳後修飾	Ubiquitinated; mediated by UHRF2 and leading to its subsequent proteasomal degradation.
細胞内局在	Nucleus.

## 画像



Anti-PCNP antibody [2707C2a] (ab70733) + immunising recombinant protein

**Predicted band size:** 19 kDa

**Observed band size:** 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors