

Product datasheet

Anti-**PLEKHG4B** antibody ab122325

画像数 2

製品の概要

製品名	Anti- <b>PLEKHG4B</b> antibody
製品の詳細	Rabbit polyclonal to <b>PLEKHG4B</b>
由来種	Rabbit
アプリケーション	<b>適用あり</b> : ICC/IF, IHC-P
種交差性	<b>交差種</b> : Human
免疫原	antigen: YDRVDEEVHR LVLTSNNRLQ QLEHLRELAS LLEGNDQQSC QKGLQLAKEN PQRTEEMVQD FRRGLSAVVS QAECREGELA RWTRSSELCE TVS, corresponding to amino acids 437-529 of Human <b>PLEKHG4B</b> (Q96PX9). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Human corpus, uterine tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122325** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment
IHC-P		1/1000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報


配列類似性	Contains 1 DH (DBL-homology) domain. Contains 1 PH domain.
-------	---

## 画像



Immunofluorescent staining of Human cell line A-431 shows positivity in nucleus and cytoplasm. Recommended concentration of ab122325 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-  
PLEKHG4B antibody (ab122325)



Immunohistochemical staining (Formalin/PFA-fixed paraffin-embedded sections) of human corpus, uterine with ab122325 at 1/1000 dilution showing strong cytoplasmic positivity in glandular cells.

Immunohistochemistry (Formalin/PFA-fixed paraffin-  
embedded sections) - Anti-PLEKHG4B antibody  
(ab122325)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors