

## Product datasheet

# Anti-POLD3 antibody ab74307

### 画像数 1

#### 製品の概要

製品名	Anti-POLD3 antibody
製品の詳細	Rabbit polyclonal to POLD3
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB <b>適用なし:</b> IP
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Rabbit, Dog, Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Bat ▲
免疫原	Synthetic peptide corresponding to a region between residues 416 and 466 of human POLD3 (NP_006582.1)
ポジティブ・コントロール	Lysate from HeLa or 293T cells

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 1.815% Tris, 1.764% Sodium citrate, 0.021% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab74307** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/10000. Predicted molecular weight: 51 kDa.

**追加情報** Is unsuitable for IP.

## ターゲット情報

機能	Required for optimal DNA polymerase delta activity.
細胞内局在	Nucleus.

## 画像



**All lanes :** Anti-POLD3 antibody (ab74307) at 0.1 µg/ml

**Lane 1 :** Whole cell lysate from HeLa cells at 50 µg

**Lane 2 :** Whole cell lysate from HeLa cells at 15 µg

**Lane 3 :** Whole cell lysate from HeLa cells at 5 µg

**Lane 4 :** Whole cell lysate from 293T cells at 50 µg

**Predicted band size:** 51 kDa

**Observed band size:** 63 kDa

[why is the actual band size different from the predicted?](#)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors