

Product datasheet

Anti-PRRT3 antibody ab122252

画像数 2

製品の概要

製品名	Anti-PRRT3 antibody
製品の詳細	Rabbit polyclonal to PRRT3
由来種	Rabbit
アプリケーション	適用あり: IHC-P, ICC/IF
種交差性	交差種: Human
免疫原	antigen sequence: PLENSEIPMI PGAHPKGSVG SEPQAFDVFP ENPRADSHRN SDVRHAPAE MPEKPVASPL GPALYGPKAA QGAQRERLPV TDDLQMAQGP SSH, corresponding to amino acids 33-125 of Human PRRT3 (UniProt Q5FWE3). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	WB: RT-4 cell and Human plasma lysates IHC-P: Human duodenum tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122252** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
----------	-----------	------

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment

ターゲット情報


細胞内局在 Membrane.

画像



Immunofluorescent staining of Human cell line U-2 OS shows positivity in plasma membrane. Recommended concentration of ab122252 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-PRRT3 antibody (ab122252)



ab122252, at a 1/250 dilution, staining PRRT3 in paraffin embedded Human duodenum tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PRRT3 antibody (ab122252)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors