

Product datasheet

Anti-PTCHD2 antibody ab113414

画像数 3

製品の概要

製品名	Anti-PTCHD2 antibody
製品の詳細	Rabbit polyclonal to PTCHD2
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Recombinant N terminal fragment of Human PTCHD2
ポジティブ・コントロール	Human Placenta, Testis and Uterus tissues.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99% PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab113414** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 µg/ml.

ターゲット情報


組織特異性	Expressed in brain and testis.
配列類似性	Belongs to the patched family. Contains 1 SSD (sterol-sensing) domain.
発生段階	Expressed in fetal brain.
細胞内局在	Membrane.

画像




ab113414 at 5 µg/ml staining PTCHD2 in Human placenta by Immunohistochemistry
Formalin-fixed, Paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTCHD2 antibody (ab113414)



ab113414 at 5 µg/ml staining PTCHD2 in Human testis by Immunohistochemistry
Formalin-fixed, Paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTCHD2 antibody (ab113414)



ab113414 at 5 µg/ml staining PTCHD2 in Human uterus by Immunohistochemistry Formalin-fixed, Paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTCHD2 antibody (ab113414)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors