


Product datasheet

Anti-Pleckstrin antibody ab17020

4 References [画像数 1](#)

製品の概要

製品名	Anti-Pleckstrin antibody
製品の詳細	Goat polyclonal to Pleckstrin
由来種	Goat
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Dog 
免疫原	Synthetic peptide: EWIKAIQMASRTGK, corresponding to C terminal amino acids 337-350 of Human Pleckstrin. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	WB: Human heart tissue and peripheral blood mononucleocyte lysates.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris buffered saline. pH 7.3
精製度	Immunogen affinity purified
特記事項(精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab17020** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリ ケーション	Abreviews	特記事項
--------------	-----------	------

WB Use a concentration of 0.05 - 1 µg/ml. Detects a band of approximately 45 kDa (predicted molecular weight: 40 kDa). Can be blocked with [Human Pleckstrin peptide \(ab23303\)](#).

ターゲット情報

機能	Major protein kinase C substrate of platelets.
配列類似性	Contains 1 DEP domain. Contains 2 PH domains.

画像



ab17020 at 0.05 µg/ml detecting Pleckstrin from Human PBMC lysate (35µg total protein per lane in RIPA buffer). Primary antibody was incubated for 1 hour. The blot was detected by chemiluminescence.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors