

Product datasheet

Anti-Poliovirus type 2 antibody [29402] ab47820

1 References

製品の概要

製品名	Anti-Poliovirus type 2 antibody [29402]
製品の詳細	Mouse monoclonal [29402] to Poliovirus type 2
由来種	Mouse
アプリケーション	適用あり: ELISA
種交差性	Reacts with Polio Virus type 2. Not yet tested in other species.
免疫原	Polio Virus type 2 (inactivated)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 0.5M Sodium chloride, 0.01M PBS, pH 7.4
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	29402
ミエローマ	x63-Ag8.653
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab47820** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

ELISA

追加情報

ELISA: Use at an assay dependent dilution.
Reacts well in ELISA both as catching antibody and as detection antibody. Weak neutralization

of Polio Virus type 2 has been observed.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性

Polio (acute poliomyelitis, infantile paralysis) is caused by the poliovirus, which may infect the nervous system and cause lasting paralysis. Poliovirus is a single stranded RNA virus and a member of the Picornavirus family. It has three strains: Poliovirus 1, 2 and 3. Replication of poliovirus RNA is accomplished by the error prone viral RNA dependent RNA polymerase. The cell receptor for poliovirus is CD155.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors