

## Product datasheet

# Anti-Prostate Specific Antigen antibody ab35355

### 1 References

#### 製品の概要

製品名	Anti-Prostate Specific Antigen antibody
製品の詳細	Sheep polyclonal to Prostate Specific Antigen
由来種	Sheep
特異性	Minimal reaction in gel diffusion with serum proteins: albumin, retinol binding protein, IgG and alpha-1-microglobulin.
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Human
免疫原	Full length native Prostrate Specific Antigen (PSA) from seminal fluid (Human).

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.2
精製度	IgG fraction
特記事項(精製)	Ammonium sulphate fractionation
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab35355** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/5000 - 1/20000.

## ターゲット情報

---

機能	Hydrolyzes semenogelin-1 thus leading to the liquefaction of the seminal coagulum.
配列類似性	Belongs to the peptidase S1 family. Kallikrein subfamily. Contains 1 peptidase S1 domain.
細胞内局在	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors