

Product datasheet

Anti-Pygopus 2 antibody [PYGAD53A] ab51367

画像数 1

製品の概要

製品名	Anti-Pygopus 2 antibody [PYGAD53A]
製品の詳細	Mouse monoclonal [PYGAD53A] to Pygopus 2
由来種	Mouse
アプリケーション	適用あり: Dot blot, ELISA, WB
種交差性	交差種: Human
免疫原	Recombinant fragment (Human)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
精製度	Protein G purified
特記事項(精製)	Filtered through a 0.22 µm membrane.
ポリ/モノ	モノクローナル
クローン名	PYGAD53A
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab51367** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 41 kDa.

ターゲット情報

機能	Involved in signal transduction through the Wnt pathway.
配列類似性	Contains 1 PHD-type zinc finger.
細胞内局在	Nucleus.

画像



Western blot - Anti-Pygopus 2 antibody
[PYGAD53A] (ab51367)

Lane 1 : Molecular weight markers

Lane 2 : Anti-Pygopus 2 antibody [PYGAD53A] (ab51367)

Lane 2 : immunizing peptide

Predicted band size: 41 kDa

Observed band size: 35,38 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors