

Product datasheet

Anti-RASA3 antibody ab31286

製品の概要

製品名 Anti-RASA3 antibody  
 製品の詳細 Goat polyclonal to RASA3

**i** This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

由来種 Goat

種交差性 交差が予測される動物種: Mouse, Rat, Dog, Human **!**

免疫原 Synthetic peptide:  
 C-QYKRDKFKKTKYGSQ  
 , corresponding to C terminal amino acids 796-810 of Human RASA3

[Run BLAST with](#) [Run BLAST with](#)

製品の特性

製品の状態 Liquid  
 保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.  
 バッファー Preservative: 0.02% Sodium Azide  
 Constituents: 0.5% BSA, Tris buffered saline. pH 7.3  
 精製度 Immunogen affinity purified  
 特記事項(精製) Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.  
 ポリ/モノ ポリクローナル  
 アイソタイプ IgG

アプリケーション

追加情報 This antibody gave a positive result in ELISA against the immunizing peptide. Antibody detection limit dilution 1:32,000.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

機能	Inhibitory regulator of the Ras-cyclic AMP pathway. Binds inositol tetrakisphosphate (IP4) with high affinity. Might be a specific IP4 receptor.
配列類似性	Contains 1 Btk-type zinc finger. Contains 2 C2 domains. Contains 1 PH domain. Contains 1 Ras-GAP domain.
細胞内局在	Cell membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors