

## Product datasheet

# Rabbit Anti-Horse IgG H&L ab102397

### 製品の概要

製品名	Rabbit Anti-Horse IgG H&L
製品の詳細	Rabbit polyclonal Secondary Antibody to Horse IgG - H&L
由来種	Rabbit
ターゲット生物種	Horse
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with horse IgG and with light chains common to other horse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species.
アプリケーション	<b>適用あり:</b> IHC-P, ICC, ELISA, WB

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab102397** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/2000.
ICC		1/200 - 1/2000.
ELISA		1/100 - 1/500. for coating; 1/1000 - 1/30000 as primary.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. colorimetric detection; 1/1000 - 1/30000 chemiluminescent detection.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors