

## Product datasheet

# Rat IgG2b, kappa monoclonal [RTK4530] (PE/Cy5®) - Isotype control ab18538

### 1 References

#### 製品の概要

製品名	Rat IgG2b, kappa monoclonal [RTK4530] (PE/Cy5®) - Isotype control
標識	PE/Cy5®. Ex: 496nm, Em: 670nm
特記事項 (標識)	The solution is free of unconjugated PE/Cy5 and unconjugated antibody.
アプリケーション	<b>適用あり:</b> Flow Cyt
免疫原	Trinitrophenol and KLH
特記事項	Isotype control has been screened on a variety of resting, activated, live, and fixed mouse and human tissues.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	モノクローナル
クローン名	RTK4530
アイソタイプ	IgG2b
軽鎖の種類	kappa

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab18538** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		

## 追加情報

FACS: Use at an assay dependent dilution for surface or intracellular staining.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors