

Product datasheet

Recombinant Cow IL17 protein ab87905

製品の概要

製品名	Recombinant Cow IL17 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Yeast

アミノ酸配列

生物種	Cow
配列	KAGVIIPQSP GCPPTEDKNF PQHVRVNLNI VNRSTNSRRP TDYHKRSTSP WTLHRNEDPE RYPSTVIWEAK CSHSGCINAE GKVDHMMNSV TIQQEILVLR RESQHCPHSF RLEKMLVAVG CTCVTPIVRH LA

特性

Our [Abpromise guarantee](#) covers the use of **ab87905** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
精製度	> 95 % SDS-PAGE. ab87905 is purified by Ion-exchange chromatography.
製品の状態	Lyophilised

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 10% Trehalose, 1X PBS
再構成	Reconstitute with sterile phosphate-buffered saline containing at least 0.1% carrier protein. Reconstituted protein is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. Avoid repeated freeze/thaw cycles.

関連情報

機能	Ligand for IL17RA and IL17RC (PubMed:17911633). The heterodimer formed by IL17A and IL17F is a ligand for the heterodimeric complex formed by IL17RA and IL17RC (PubMed:18684971). Involved in inducing stromal cells to produce proinflammatory and hematopoietic cytokines (PubMed:8676080).
組織特異性	Restricted to activated memory T-cells.
配列類似性	Belongs to the IL-17 family.
翻訳後修飾	Found both in glycosylated and nonglycosylated forms.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors