

Product datasheet

Recombinant HIV1 tat protein ab83353

1 References

製品の概要

製品名	Recombinant HIV1 tat protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
領域	1 to 101

特性

Our [Abpromise guarantee](#) covers the use of **ab83353** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生理活性	ab83353 reacts with anti-Tat antibodies from human, monkey, rabbit and mouse serum.
アプリケーション	Blocking Western blot SDS-PAGE
精製度	> 90 % SDS-PAGE.
製品の状態	Lyophilised
備考	ab83353 reacts with anti-Tat antibodies from human, monkey, rabbit and mouse serum.

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -80°C. Preservative: None Constituents: 0.1% Glycerol
再構成	The protein should be reconstituted in apyrogenic sterile water or PBS buffer. The reconstituted solution has to be used immediately.

関連情報

関連性	HIV1 tat (Transactivator of Transcription) protein is a pleiotropic factor that induces a broad range of biological effects in numerous cell types. At the HIV promoter, tat is a powerful transactivator of gene expression which acts by inducing chromatin remodeling and by recruiting elongation-competent transcriptional complexes on to the viral LTR.
細胞内局在	Nuclear. Upon stimulation with EDN1, it is exported from the nucleus to the perinuclear region.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors