

Product datasheet

Recombinant Human Annexin A10 protein ab123206

画像数 1

製品の概要

| | |
|--------|---------------------------------------|
| 製品名 | Recombinant Human Annexin A10 protein |
| タンパク質長 | Full length protein |

製品の詳細

| | |
|-----------|------------------------|
| 由来 | Recombinant |
| 由来 | Escherichia coli |
| アミノ酸配列 | |
| アクセッション番号 | Q9UJ72 |
| 生物種 | Human |

| | |
|----|---|
| 配列 | MGSSHHHHHH SGLVPRGSH MGSFMFCGDY VQGTIFPAPN FNPIMDAQML GGALQGFDCD KDMLINILTQ RCNAQRMMIA EAYQSMYGRD LIGDMREQLS DHFKDVMAGL MYPPPLYDAH ELWHAMKGVG TDENCLIEIL ASRTNGEIFQ MREAYCLQYS NNLQEDIYSE TSGHFRDTLM NLVQGTREEG YTDPAQAQD AMVLWEACQQ KTGEHKTMLQ MILCNKSYQQ LRLVFQEFQN ISGQDMVDAI NECYDGYFQE LLVAIVLCVR DKPAYFAYRL YSAIHDFGFH NKTVIRILIA RSEIDLLTIR KRYKERYGKS LFHDIRNFAS GHYKKALLAI CAGDAEDY |
|----|---|

| | |
|-----|-----------------------|
| 分子量 | 40 kDa including tags |
| 領域 | 1 to 324 |
| タグ | His tag N-Terminus |

特性

Our [Abpromise guarantee](#) covers the use of **ab123206** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|-------------------|
| アプリケーション | Mass Spectrometry |
| | SDS-PAGE |

| | |
|-------|--|
| 質量分析 | MALDI-TOF |
| 精製度 | > 95 % SDS-PAGE. ab123206 is purified by using conventional chromatography. |
| 製品の状態 | Liquid |


前処理および保存

| | |
|------------|--|
| 保存方法および安定性 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol, 0.58% Sodium chloride |
|------------|--|

関連情報

| | |
|-------|---|
| 配列類似性 | Belongs to the annexin family. Contains 4 annexin repeats. |
|-------|---|

画像



15% SDS-PAGE analysis of ab123206 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors