

## Product datasheet

# Recombinant Human COUP TF1 protein ab112360

### 画像数 1

#### 製品の概要

製品名	Recombinant Human COUP TF1 protein
タンパク質長	Protein fragment

#### 製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	<a href="#">P10589</a>
生物種	Human
配列	CRLKKCLKVGMRRREAVQRGRMPPTQPNPGQYALTNGDPLN GHCYLSGYIS LLLRAEPYPTSRYGSSQCMQPNNIMGIENICELARLLFSA VEWARNIPFF
分子量	37 kDa including tags
領域	141 to 240

#### 特性

Our [Abpromise guarantee](#) covers the use of **ab112360** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Western blot ELISA
製品の状態	Liquid
備考	Best use within three months from the date of receipt of this protein.

#### 前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------	---

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## 関連情報

機能	Coup (chicken ovalbumin upstream promoter) transcription factor binds to the ovalbumin promoter and, in conjunction with another protein (S300-II) stimulates initiation of transcription. Binds to both direct repeats and palindromes of the 5'-AGGTCA-3' motif. Represses transcriptional activity of LHCG.
配列類似性	Belongs to the nuclear hormone receptor family. NR2 subfamily. Contains 1 nuclear receptor DNA-binding domain.
細胞内局在	Nucleus.

## 画像



ab112360 on a 12.5% SDS-PAGE Stained with Coomassie Blue.

SDS-PAGE - Recombinant Human COUP TF1  
protein (ab112360)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors