

Product datasheet

Recombinant Human EPB41L1 protein ab117026

画像数 1

製品の概要

製品名	Recombinant Human EPB41L1 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	Q9H4G0-2
生物種	Human
配列	MEEKDYSEADGLSERTTPSKAQKSPQKIAKKYKSAICRVTLDDASEYECE VEKHGRGQVLFDLVCEHLNLLLEKDYFGLTFCDADSQKN
分子量	35 kDa including tags
領域	1 to 88

特性

Our [Abpromise guarantee](#) covers the use of **ab117026** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE ELISA
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------	---

Constituents: 0.79% Tris HCl, 0.3% Glutathione

関連情報

機能	May function to confer stability and plasticity to neuronal membrane via multiple interactions, including the spectrin-actin-based cytoskeleton, integral membrane channels and membrane-associated guanylate kinases.
組織特異性	Highest expression in brain, lower in heart, kidney, pancreas, placenta, lung and skeletal muscle.
配列類似性	Contains 1 FERM domain.
細胞内局在	Cytoplasm > cytoskeleton.

画像



12.5% SDS-PAGE showing ab117026 at approximately 35.31 kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors