

Product datasheet

Recombinant Human Glycerol 3 Phosphate Dehydrogenase protein ab93724

画像数 1

製品の概要

製品名	Recombinant Human Glycerol 3 Phosphate Dehydrogenase protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
生物種	Human

配列	MASKKVCIVG SGNWGSIAIK IVGGNAAQLA
	QFDPRTVMWV FEEDIGGKKL TEIINTQHEN
	VKYLPGHKLP PNVVAVPDVV QAAEDADILI
	FVVPHQFIGK ICDQLKGHLK ANATGISLIK
	GVDEGPNGLK LISEVIGERL GIPMSVLMGA
	NIASEVADEK FCETTIGCKD PAQGQLLKEL
	MQTPNFRITV VQEVDTVEIC GALKNVVAVG
	AGFCDGLGFG DNTKAAVIRL GLMEMIAFAK
	LFCSGPVSSA TFLESCGVAD LITTCYGGRN
	RKVAEAFART GKSIEQLEKE LLNGQKLQGP
	ETARELYSIL QHKGLVDKFP LFMAVYKVCY
	EGQPVGFEIH CLQNHPEHM

特性

Our [Abpromise guarantee](#) covers the use of **ab93724** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
精製度	> 90 % SDS-PAGE. ab93724 is purified using conventional chromatography techniques.
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.316% Tris HCl, 10% Glycerol

関連情報


配列類似性

Belongs to the NAD-dependent glycerol-3-phosphate dehydrogenase family.

細胞内局在

Cytoplasm.

画像



15% SDS-PAGE showing ab93724 at approximately 37.5kDa (3µg).

SDS-PAGE - Recombinant Human Glycerol 3
Phosphate Dehydrogenase protein (ab93724)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors