

## Product datasheet

# Recombinant Human HDHD3 protein ab125485

### 製品の概要

製品名	Recombinant Human HDHD3 protein
タンパク質長	Full length protein

### 製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列  
アクセッション番号

[Q9BSH5](#)

生物種 Human

配列

```
MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD
KDRWGSMahr LQIRLLTWDV KDTLLRLRHP
LGEAYATKAR AHGLEVEPSA LEQGFRQAYR
AQSHSFPNYG LSHGLTSRQW WLDVVLQTFH
LAGVQDAQAV APIAEQLYKD FSHPCTWQVL
DGAEDTLREC RTRGLRLAVI SNFDRRLEGI
LEGLGLREHF DFVLTSEAAG WPKPDPRIFF
EALRLAHMEP VVAHVGDNY LCDYQGPRAV
GMHSFLVVG P QALDPVVRDS VPKEHILPSL
AHL L PALDCL EGSTPGL
```

分子量 32 kDa including tags

領域 1 to 251

タグ His tag N-Terminus

配列の追加情報 HDHD3 Human Recombinant fused with a36 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 287 amino acids (1-251 a.a.).

### 特性

Our [Abpromise guarantee](#) covers the use of **ab125485** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

精製度 > 95 % SDS-PAGE.

ab125485 is purified by proprietary chromatographic techniques.

製品の状態	Liquid
備考	For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

---

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Store under desiccating conditions. pH: 8.00 Constituents: 0.32% Tris HCl, 10% Glycerol
------------	---

---

#### 関連情報

配列類似性	Belongs to the HAD-like hydrolase superfamily.
-------	--

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

---

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors