

## Product datasheet

# Recombinant Human HIGD1A protein ab125988

### 製品の概要

---

製品名	Recombinant Human HIGD1A protein
タンパク質長	Full length protein

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

アクセッション番号	<a href="#">Q9Y241</a>
生物種	Human
配列	MSTDTGVSLPSYEEDQGSKLIRKAKEAPFVPGIAGFAAIVAYGLYKLS RGNTKMSIHLIHMVAAQGFVVGAMTVGMGYSMYREFWAKPKP
領域	1 to 93

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab125988** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilised

### 前処理および保存

---

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
再構成	Reconstitute with water to desired concentration.

### 関連情報

---

配列類似性	Contains 1 HIG1 domain.
細胞内局在	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors