

Product datasheet

Recombinant Human LBP protein ab114808

画像数 1

製品の概要

製品名	Recombinant Human LBP protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	<a href="#">P18428</a>
生物種	Human
配列	<p>ANPGLVARITDKGLQYAAQEGLLALQSELLRITLPDFTGD            LRIPHVGRGR            YEFHSLNIHSCCELLHSALRPVPGQLSLSISDSSIRVQGR            WKVRKSFFKL            QGSFDVSVKGISISVDLLLGSSESSGRPTVTASSCSSDIAD            VEVDMSGDLG            WLLNLFHNQIESKFQKVLESRICEMIQSVSSDLQPYLQT            LPVTTEIDSF            ADIDYSLVEAPRATAQMLEVMFKGEIFHRNHRSPVTLAA            VMSLPEEHNK            MVYFAISDYVFNTASLVYHEEGLNFSITDDMIPPDSNIR            LTTKSFRPFV            PRLARLYPNMNEELQGSVPSAPLLNFSPGNLSVDPYMEID            AFVLLPSSSK            EPVFRLSVATNVSATLTFNTSKITGFLKPGKVKVELKESK            VGLFNAELLE            ALLNYYILNTFYPKFNDKLAEGFPLPLLKRQYLDLGLQI            HKDFPSRRDH PP</p>
分子量	76 kDa including tags
領域	26 to 477

特性

Our [Abpromise guarantee](#) covers the use of **ab114808** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション**

Western blot

SDS-PAGE

ELISA

**製品の状態**

Liquid

**備考**

Protein concentration is above or equal to 0.05 mg/ml.

Best used within three months from the date of receipt.

**前処理および保存**

**保存方法および安定性**

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

**関連情報**

**機能**

Binds to the lipid A moiety of bacterial lipopolysaccharides (LPS), a glycolipid present in the outer membrane of all Gram-negative bacteria. The LBP/LPS complex seems to interact with the CD14 receptor.


**配列類似性**

Belongs to the BPI/LBP/Plunc superfamily. BPI/LBP family.

**細胞内局在**

Secreted.

**画像**



ab114808 on 12.5% SDS-PAGE Stained with Coomassie Blue.

SDS-PAGE - Recombinant Human LBP protein  
(ab114808)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors