

Product datasheet

Recombinant Human NME3 protein ab103051

画像数 1

製品の概要

製品名	Recombinant Human NME3 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
アクセッション番号	Q13232
生物種	Human

配列	<p>MGSSHHHHHSSGLVPRGSHMERTFLAVKPDGVQRRLV GEIVRRFERKGF KLVALKLVQASEELLREHYAELRERPFYGRLVKYMASGPV VAMVWQGLDV VRTSRALIGATNPADAPPGTIRGDFCIEVGKNLIHGSDSV ESARREIALW FRADELLCWEDSAGHWLYE</p>
----	--

分子量	19 kDa including tags
領域	22 to 169
タグ	His tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab103051** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	<p>SDS-PAGE</p> <p>Mass Spectrometry</p>
質量分析	MALDI-TOF
精製度	<p>> 95 % SDS-PAGE.</p> <p>ab103051 was purified using conventional chromatography.</p>
製品の状態	Liquid

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.0308% DTT, 0.316% Tris HCl, 50% Glycerol, 0.58% Sodium chloride

関連情報

機能

Major role in the synthesis of nucleoside triphosphates other than ATP. The ATP gamma phosphate is transferred to the NDP beta phosphate via a ping-pong mechanism, using a phosphorylated active-site intermediate. Probably has a role in normal hematopoiesis by inhibition of granulocyte differentiation and induction of apoptosis.


配列類似性

Belongs to the NDK family.

発生段階

Preferentially expressed at early stages of myeloid differentiation of highly purified CD34+ cells.

画像



15% SDS-PAGE analysis of 3µg ab103051

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors