

Product datasheet

Recombinant Human PC6 protein ab114944

画像数 1

製品の概要

製品名	Recombinant Human PC6 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	Q92824
生物種	Human
配列	TEGYFMEDGRCVQSCSISYYFDHSSSENGYKSCCKCDISCL TCNGPGFKNC TSCPSGYLLDLGMCQMGAICKDATEESWAEGGFCMLVKKN NLCQRKVLQQ LCCKTCTFQG
分子量	38 kDa including tags
領域	804 to 913

特性

Our [Abpromise guarantee](#) covers the use of **ab114944** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best used within three months from the date of receipt.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------	---


pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能	Likely to represent a widespread endoprotease activity within the constitutive and regulated secretory pathway. Capable of cleavage at the RX(K/R)R consensus motif.
組織特異性	Expressed in T-lymphocytes.
配列類似性	Belongs to the peptidase S8 family. Contains 1 homo B/P domain. Contains 1 PLAC domain.
ドメイン	The propeptide domain acts as an intramolecular chaperone assisting the folding of the zymogen within the endoplasmic reticulum.
細胞内局在	Secreted.

画像



SDS-PAGE - Recombinant Human PC6 protein
(ab114944)

ab114944 analyzed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors