

Product datasheet

Recombinant Human PFKFB1 protein ab114752

画像数 1

製品の概要

製品名	Recombinant Human PFKFB1 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	P16118
生物種	Human
配列	<p>MSPEMGELTQTRLQKIWIIPHSSGSSRLQRRRGSSIPQFTN SPTMVIMVGL PARGKTYISTKLTRYLNWIGTPTKVFNLGQYRREAVSYKN YEFFLPDNME ALQIRKQCALAALKDVHNYLSHEEGHVAVFDATNTTRERR SLILQFAKEH GYKVFFIESICNDPGIIAENIRQVKGSPDYIDCDREKVL EDFLKRIECY EVNYQLDEELDSHLSYIKIFDVGTRYMVNRVQDHIQSRT VYYLMNIHVT PRSIYLCRHGESELNIRGRIGGDSGLSVRGKQYAYALANF IQSQGISSLK VWTSRMKRTIQTAEALGVPYEQWKALNEIDAGVCEEMTYE EIQEHYPEEF ALRDQDKYRYRYPKGESYEDLVQRLEPVIMELERQENVLV ICHQAVMRCL LAYFLDKSSDELPLYLKCPLHTVLKLTVPVAYGCKVESIYLN VEAVNTHREK PENVDITREPEEALDTPAHY</p>
分子量	78 kDa including tags
領域	1 to 471

特性

Our [Abpromise guarantee](#) covers the use of **ab114752** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション

ELISA
SDS-PAGE
Western blot

製品の状態

Liquid

備考

Protein concentration is above or equal to 0.05 µg/ul.
Best use within three months from the date of receipt of this protein. EC=3.1.3.46

前処理および保存

保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.3% Glutathione, 0.79% Tris HCl


関連情報

機能 Synthesis and degradation of fructose 2,6-bisphosphate.

組織特異性 Liver.

配列類似性 In the C-terminal section; belongs to the phosphoglycerate mutase family.

画像



12.5% SDS-PAGE Stained with Coomassie Blue showing ab114752 at approximately 77.88kDa.

SDS-PAGE - Recombinant Human PFKFB1 protein (ab114752)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors