

Product datasheet

Recombinant Human RNPEP protein ab114879

画像数 1

製品の概要

製品名	Recombinant Human RNPEP protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	Q9H4A4
生物種	Human
配列	GNVKKLGDITYPSISNARNAELRLRWGQIVLKNDHQEDFWKVKEFLHNQGK QKYTLPLYHAMMGSEVAQTLAKETFASTASQLHSNVVNYVQQIVAPKGS
分子量	37 kDa including tags
領域	551 to 650

特性

Our [Abpromise guarantee](#) covers the use of **ab114879** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------	---

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能	Exopeptidase which selectively removes arginine and/or lysine residues from the N-terminus of several peptide substrates including Arg(0)-Leu-enkephalin, Arg(0)-Met-enkephalin and Arg(-1)-Lys(0)-somatostatin-14. Can hydrolyze leukotriene A4 (LTA-4) into leukotriene B4 (LTB-4).
配列類似性	Belongs to the peptidase M1 family.
細胞内局在	Secreted.

画像



SDS-PAGE - RNPEP protein (ab114879)

12.5% SDS-PAGE showing ab114879 at approximately 36.63kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors