

Product datasheet

Recombinant Human Rab2 protein ab116399

画像数 1

製品の概要

| | |
|--------|--------------------------------|
| 製品名 | Recombinant Human Rab2 protein |
| タンパク質長 | Full length protein |

製品の詳細

| | |
|-----------|------------------------|
| 由来 | Recombinant |
| 由来 | Escherichia coli |
| アミノ酸配列 | |
| アクセッション番号 | P61019 |
| 生物種 | Human |
| 分子量 | 53 kDa including tags |
| 領域 | 1 to 212 |

特性

Our [Abpromise guarantee](#) covers the use of **ab116399** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|--|
| アプリケーション | Western blot ELISA SDS-PAGE Mass Spectrometry |
| 質量分析 | MALDI-TOF-TOF |
| 精製度 | > 95 % SDS-PAGE. |
| 製品の状態 | Liquid |


前処理および保存

| | |
|------------|---|
| 保存方法および安定性 | Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.40 Constituent: 99% PBS |
|------------|---|

関連情報

| | |
|-------|--|
| 機能 | Required for protein transport from the endoplasmic reticulum to the Golgi complex. |
| 配列類似性 | Belongs to the small GTPase superfamily. Rab family. |
| 細胞内局在 | Endoplasmic reticulum-Golgi intermediate compartment membrane. Melanosome. Endoplasmic reticulum membrane. Golgi apparatus membrane. Identified by mass spectrometry in melanosome fractions from stage I to stage IV. |

画像



10% SDS-PAGE stained with Coomassie Blue (Lane 1),
immunoblotting with anti-6xHis (Lane 2)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors