

## Product datasheet

# Recombinant Human Ran (mutated T24 N) protein ab90206

### 製品の概要

---

製品名	Recombinant Human Ran (mutated T24 N) protein
タンパク質長	Full length protein

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

生物種	Human
修飾	mutated T24N

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab90206** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
精製度	> 95 % SDS-PAGE. Protein preparation is 18% GDP-loaded, measured by HPLC.

製品の状態	Liquid
-------	--------

備考	Ran T24N mutation results in no or little binding to GTP and is known to act as a negative inhibitor of nuclear import.
----	---

### 前処理および保存

---

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Constituents: 0.1001% DTE (1,4-Dithioerythritol), 1.125% Potassium chloride, 0.42% Potassium phosphate, 0.0475% Magnesium chloride, 10% Glycerol
------------	---

### 関連情報

<b>機能</b>	GTP-binding protein involved in nucleocytoplasmic transport. Required for the import of protein into the nucleus and also for RNA export. Involved in chromatin condensation and control of cell cycle (By similarity). The complex with BIRC5/ survivin plays a role in mitotic spindle formation by serving as a physical scaffold to help deliver the RAN effector molecule TPX2 to microtubules. Enhances AR-mediated transactivation. Transactivation decreases as the poly-Gln length within AR increases.
<b>組織特異性</b>	Expressed in a variety of tissues.
<b>配列類似性</b>	Belongs to the small GTPase superfamily. Ran family.
<b>翻訳後修飾</b>	The N-terminus is blocked.
<b>細胞内局在</b>	Nucleus. Cytoplasm. Melanosome. Becomes dispersed throughout the cytoplasm during mitosis. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors