

## Product datasheet

# Recombinant Mouse IL17E protein ab119730

### 製品の概要

---

製品名	Recombinant Mouse IL17E protein
タンパク質長	Full length protein

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

アクセッション番号 [Q8VHH8](#)

生物種 Mouse

配列  
 VSLRIQEGCSHL PSCCPSKEQPPEEWLKWSSASVSPPEPLSHTTHAESC  
 RASKDGPLNSRAISPWSYELDRDLNRVPQDLYHARCLCPHCVSLQTGSHM  
 DPLGNSVPLYHNQTVFYRRPCHGEEGTRRYCLERRLYRVSLACVCVRPRV  
 MA

分子量 35 kDa

領域 17 to 169

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab119730** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	HPLC Functional Studies
精製度	> 95 % SDS-PAGE. Purity by SDS-PAGE and HPLC: = or > 95% Endotoxin level is < 0.01 ng per µg of ab119730.
製品の状態	Lyophilised
備考	The lyophilized ab119730 is best-stored desiccated at -20 °C. Upon reconstitution the product should be aliquoted and stored at -20 °C. It is recommended to add a carrier protein (0.1% HSA or BSA) for long term storage.

### 前処理および保存

---

<b>保存方法および安定性</b>	Shipped at 4°C. Store at -20°C.
<b>再構成</b>	Centrifuge vial before opening. When reconstituting the product, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution. It is recommended to reconstitute the lyophilized product with sterile 5 mM Acetic acid to a concentration of 0.1 mg/ml, which can be further diluted into other aqueous solutions.
<b>関連情報</b>	
<b>機能</b>	Induces activation of NF-kappa-B and stimulates production of the proinflammatory chemokine IL-8. Proinflammatory cytokine favoring Th2-type immune responses.
<b>組織特異性</b>	Expressed at low levels in several tissues, including brain, kidney, lung, prostate, testis, spinal cord, adrenal gland, and trachea.
<b>配列類似性</b>	Belongs to the IL-17 family.
<b>細胞内局在</b>	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors