

Recombinant Mouse IL-28A protein ab50172

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant Mouse IL-28A protein |
| 精製度 | > 95 % SDS-PAGE. ab50172 purity was assessed also by HPLC Endotoxin level is less than 0.1 ng per µg (1EU/µg). |
| エンドトキシン・レベル | < 0.100 Eu/µg |
| 発現系 | Escherichia coli |
| タンパク質長 | Protein fragment |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Mouse |
| 配列 | MDPVPRATRL PVEAKDCHIA QFKSLSPKEL QAFKKAKDAI EKRLLEKDMR CSSHLISRAW DLKQLQVQER PKALQAEVAL TLKVWENMTD SALATILGQP LHTLSHIHSQ LQTCTQLQAT AEPKPPSRRL SRWLHRLQEA QSKETPGCLE DSVTSNLFRL LTRDLKCVAS GDQCV |
| 予測される分子量 | 20 kDa |
| 領域 | 19 to 193 |

特性

Our **Abpromise guarantee** covers the use of **ab50172** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|--------------------------------|
| アプリケーション | SDS-PAGE Functional Studies |
|----------|--------------------------------|

| | |
|-------|-------------|
| 製品の状態 | Lyophilized |
|-------|-------------|

前処理および保存

| | |
|------------|--|
| 保存方法および安定性 | Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term. |
| 再構成 | Reconstitute to 1mg/ml in distilled water. |

関連情報

| | |
|-------|--|
| 機能 | Cytokine with immunomodulatory activity. Up-regulates MHC class I antigen expression. Displays potent antiviral activity. Also displays antitumor activity. Ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway. |
| 配列類似性 | Belongs to the IL-28/IL-29 family. |
| 細胞内局在 | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors