

Product datasheet

Anti-Resistin antibody ab119501

1 References 画像数 1

製品の概要

製品名	Anti-Resistin antibody
製品の詳細	Rabbit polyclonal to Resistin
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Mouse
免疫原	Synthetic peptide corresponding to C terminal amino acids of Mouse Resistin.
ポジティブ・コントロール	Recombinant Mouse Resistin protein.

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservatives: 0.025% Thimerosal (merthiolate), 0.025% Sodium azide Constituents: 2.5% BSA, 0.45% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab119501** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリ ケーション	Abreviews	特記事項
--------------	-----------	------


WB Use a concentration of 1 µg/ml. Predicted molecular weight: 12 kDa. The detection limit for Resistin is approximately 2.5ng/lane under non-reducing and reducing conditions.

IHC-P Use a concentration of 1 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

ターゲット情報

機能	Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.
組織特異性	Expressed only in fatty tissues.
配列類似性	Belongs to the resistin/FIZZ family.
細胞内局在	Secreted.

画像



All lanes : Anti-Resistin antibody (ab119501)
at 1 µg/ml

Lane 1 : Recombinant Mouse Resistin protein
at 0.01 µg

Lane 2 : Recombinant Mouse Resistin protein
at 0.005 µg

Predicted band size: 12 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors