

Product datasheet

Anti-Respiratory Syncytial Virus antibody (Biotin) ab19986

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Respiratory Syncytial Virus antibody (Biotin) |
| 製品の詳細 | Goat polyclonal to Respiratory Syncytial Virus (Biotin) |
| 由来種 | Goat |
| 標識 | Biotin |
| アプリケーション | 適用あり: Neutralising, IHC-P |
| 免疫原 | Human isolate of RSV (confirmed). |
| 特記事項 | Covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution. |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | pH: 7.20 Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS |
| 精製度 | IgG fraction |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab19986** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|--------------|-----------|-------------------------------------|
| Neutralising | | Use at an assay dependent dilution. |
| IHC-P | | Use at an assay dependent dilution. |

ターゲット情報

関連性 Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, enveloped RNA virus. The virion is variable in shape and size with average diameter of between 120 and 300 nm. The 63 kD RSV fusion protein of the RSS 2 strain (subtype A) directs fusion of viral and cellular membranes, results in viral penetration, and can direct fusion of infected cells with adjoining cells, resulting in the formation of syncytia or multi nucleated giant cells.

細胞内局在 Virion. Host cytoplasm

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors