

## Product datasheet

# Anti-SCRN3 antibody ab117709

### 画像数 2

#### 製品の概要

製品名	Anti-SCRN3 antibody
製品の詳細	Rabbit polyclonal to SCRN3
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA, IHC-P, ICC
種交差性	交差種: Human 交差が予測される動物種: Mouse
免疫原	Synthetic peptide conjugated to KLH, corresponding to an internal region of Human SCRN3
ポジティブ・コントロール	Human breast or tonsil tissue.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab117709** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 48 kDa.
ELISA		Use a concentration of 1 - 2 µg/ml.


アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 µg/ml.
ICC		Use a concentration of 20 µg/ml.

## ターゲット情報

### 配列類似性


Belongs to the peptidase C69 family. Secernin subfamily.

## 画像



Human Breast tissue, Formalin-Fixed, Paraffin-Embedded, stained with ab117709 at 5µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SCRIN3 antibody (ab117709)



Human tonsil tissue, Formalin-Fixed, Paraffin-Embedded, stained with ab117709 at 5µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SCRIN3 antibody (ab117709)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors