

Product datasheet

Anti-SF3B4 antibody ab11803

★★★★☆ 1 Abreviews 4 References 画像数 1

製品の概要

製品名	Anti-SF3B4 antibody
製品の詳細	Goat polyclonal to SF3B4
由来種	Goat
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Xenopus laevis, Drosophila melanogaster, Zebrafish
免疫原	Synthetic peptide: AAGPISERNQDAT, corresponding to amino acids 2-14 of Human SF3B4. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Hela lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: Tris buffered saline, 0.5% BSA
精製度	Immunogen affinity purified
特記事項(精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab11803** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB ★★★★☆ Use a concentration of 0.5 - 2 µg/ml. Detects a band of approximately 50 kDa (predicted molecular weight: 49 kDa). Can be blocked with [Human SF3B4 peptide \(ab23211\)](#).

ターゲット情報

画像



Anti-SF3B4 antibody (ab11803)

Predicted band size: 49 kDa

Ab11803 staining (0.1 µg/ml) of Hela lysate (RIPA buffer, 35 µg total protein per lane).

Primary incubated for 1 hour. Detected by western blot using chemiluminescence.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors